



Anti-HA(B/Texas/6/2011)

Product Name	Anti-HA (B/Texas/6/2011)
Catalog No	IT-003-B12M3
Source	Mouse monoclonal antibody, clone 3C3
Isotype	IgG2b
Formulation	Each vial contains 100 µg purified IgG in 100 µl PBS with 40% glycerol.
Immunogen	DNA vaccine expressing HA1 (B/Texas/6/2011) protein (IT-003-B12p)
Applications	WB, ELISA, IP, IF etc.
Epitope	Linear
Specificity	Pan reactive to HA(B). No cross-reactivity to HA(A) subtype as tested.
Storage	Stable for 6 months at 4°C from the date of shipment. For long term storage, aliquot and freeze at -80°C. Non-hazardous. No MSDS required.
Use Limitation	For research use only, not for use in diagnostic procedures.

